



Curlin Pump Administration Procedure (Intermittent or Continuous Mode)

Supplies:

- Medication Bag
- Curlin Tubing; change every 24 hours
- Sodium Chloride Syringes (2)
- Alcohol Wipes
- 2 C-Sized Batteries
- Gloves
- Curlin IV Pump

Procedure:

1. Gather equipment. Clean work area. Wash hands. Wear gloves.
2. Gather (1) Sodium Chloride flush as instructed by your nurse.
3. Remove the Curlin tubing from package.
4. Twist and remove breakaway tab from tubing (yellow circle).
5. Remove the protective tab from the entry port of the medication bag.
6. Remove the cover from the spike on tubing and insert the spike into the medication bag using a pushing, twisting motion. **DO NOT TOUCH SPIKE.**
7. Squeeze the yellow piece between the thumb and forefinger to allow fluid to flow through tubing (priming) until all the air is removed from the tubing set. Release the yellow piece to stop flow.
8. Insert 2 C-sized batteries into the Curlin pump. Change the batteries daily with bag change.
9. Insert tubing set into the Curlin pump. Open top door and insert the blue tubing guide on the right side of pump. The blue arrow on the pump indicates the direction and location for the blue tubing guide.
10. Insert the yellow colored tubing guide into the cavity on the left side of the pump so the tubing lies in the middle, then push the tubing slightly into the air detector slot. The yellow arrow on the pump's surface helps to locate the slot.
11. Close door.
12. Turn the pump on by pressing the orange on/off button.

13. The pump will run through a series of tests with three beeps.
14. Press “Yes” to program.
15. Press “Yes” to repeat Rx (all pumps come preprogrammed).
16. Clean the injection cap on the end of your catheter with an alcohol wipe vigorously for 15 seconds.
17. Flush your catheter with Sodium Chloride as instructed by your nurse.
18. Clean the injection cap on the end of your catheter with a new alcohol wipe vigorously for 15 seconds.
19. Attach the Curlin tubing to the end of your injection cap.
20. Open all clamps; press “Run” to start.
21. If the green light is flashing, the pump is running.

At Completion of Infusion:

1. When infusion is completed, the display will say, “Infusion Complete.” Press “Pause,” then “Off”.
2. Close the clamp on tubing.
3. Disconnect the tubing from your injection cap.
4. Clean the injection cap on the end of your catheter with an alcohol wipe vigorously for 15 seconds then flush your catheter with Sodium Chloride as instructed by your nurse.
5. Clean the injection cap on the end of your catheter again with a new alcohol wipe vigorously for 15 seconds then flush your catheter with Heparin as instructed by your nurse.
6. Open the top door and pull the tubing out of pump. Discard bag and tubing in your household trash.